


Viper Plus GT Wheelchair

Deluxe High-Strength, Lightweight, Dual-Axle



HCPCS
K0004

New “Built-in” Seat Rail Extension System and Extendable Upholstery



- Built-in seat rail extensions and extendable upholstery easily adjusts seat depth from 16" to 18". (Patent pending) (Figure A - C)
- All-Aluminum frame.
- Dual Axle allows for multiple seat-to-floor height positions.
- Adjustable height, flip-back, removable, padded arms.
- Weighs under 32 lbs. in basic configuration (excluding front rigging).
- Over-center cross brace eliminates need for seat guides. Works with seating and positioning accessories.
- Adjustable angle back easily adjusts from 5° – 20°.
- Adjustable height back adjusts from 17" – 19" in 1" increments, adaptable for individual patient needs.
- Durable, heavy-gauge flame-retardant nylon upholstery with hook and loop fastener flaps.
- Quick release, 24" composite Mag-style wheels are lightweight and maintenance free.
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability.
- Adjustable angle caster forks have three height adjustments and angle adjustability.
- Tool-free adjustable swing-away footrests also feature easy-to-access release handles.
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight.
- Adjustable anti-tipper with flip-up wheels. Never needs to be removed.
- Comes with push-to-lock wheel locks.
- Patent Pending.

Item #	Seat Width	Armrest	Product Front Riggings	Product Seat Extension	UOM
PLA416FBFAARAD-SF	16"	Flip Back, Detachable, Adj. Height Full Arm	Swing-away Footrests	Built-in	1/cs
PLA416FBFAARAD-ELR	16"	Flip Back, Detachable, Adj. Height Full Arm	Elevating Leg rests	Built-in	1/cs
PLA416FBDAARAD-SF	16"	Flip Back, Detachable, Adj. Height Desk Arm	Swing-away Footrests	Built-in	1/cs
PLA416FBDAARAD-ELR	16"	Flip Back, Detachable, Adj. Height Desk Arm	Elevating Leg rests	Built-in	1/cs
PLA418FBDAARAD-SF	18"	Flip-Back, Detachable, Adj. Height Desk Arm	Swing-away Footrests	Built-in	1/cs
PLA418FBDAARAD-ELR	18"	Flip Back, Detachable, Adj. Height Desk Arm	Elevating Leg rests	Built-in	1/cs
PLA418FBFAARAD-SF	18"	Flip Back, Detachable, Adj. Height Full Arm	Swing-away Footrests	Built-in	1/cs
PLA418FBFAARAD-ELR	18"	Flip Back, Detachable, Adj. Height Full Arm	Elevating Leg rests	Built-in	1/cs
PLA420FBDAARAD-SF	20"	Flip Back, Detachable, Adj. Height Desk Arm	Swing-away Footrests	Built-in	1/cs
PLA420FBDAARAD-ELR	20"	Flip Back, Detachable, Adj. Height Desk Arm	Elevating Leg rests	Built-in	1/cs
PLA420FBFAARAD-SF	20"	Flip Back, Detachable, Adj. Height Full Arm	Swing-away Footrests	Built-in	1/cs
PLA420FBFAARAD-ELR	20"	Flip Back, Detachable, Adj. Height Full Arm	Elevating Leg rests	Built-in	1/cs
PLA422FBDAAR-SF	22"	Flip Back, Detachable, Adj. Height Desk Arm	Swing-away Footrests	N.A.	1/cs
PLA422FBDAAR-ELR	22"	Flip Back, Detachable, Adj. Height Desk Arm	Elevating Leg rests	N.A.	1/cs
PLA422FBFAAR-SF	22"	Flip Back, Detachable, Adj. Height Full Arm	Swing-away Footrests	N.A.	1/cs
PLA422FBFAAR-ELR	22"	Flip Back, Detachable, Adj. Height Full Arm	Elevating Legrests	N.A.	1/cs

specifications

Width Open with 16" Seat	24"	Product Weight with 18" Seat	44 lbs.
Width Open with 18" Seat	26"	Product Weight with 20" Seat	45 lbs.
Width Open with 20" Seat	28"	Product Weight with 22" Seat	46 lbs.
Width Open with 22" Seat	30"	Carton Shipping Weight with 16" Seat	51 lbs.
Seat Width Closed	12.5"	Carton Shipping Weight with 18" Seat	52 lbs.
Height	36"	Carton Shipping Weight with 20" Seat	53 lbs.
Seat-To-Floor Height	17.5" - 19.5"	Carton Shipping Weight with 22" Seat	54 lbs.
Back Height	17" - 19"	Weight Limit with 16", 18" and 20" Seat	300 lbs.
Seat Depth (16"-20" Wide Chair) (22" Wide Chair)	16" - 18" 18"	Weight Limit with 22" Seat	350 lbs.
Length with Riggings	42"	Warranty	Limited Lifetime
Product Weight with 16" Seat	43 lbs.		

accessories

Universal Forward Stabilizer

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

Anti-Tippers

STDS807	Anti-Tippers With Wheels	1/pr
---------	--------------------------	------

Seat belts

STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric	1/ea

Telescoping I.V. Pole Attachment

STDS820	Telescoping IV Pole Attachment	1/ea
---------	--------------------------------	------

Universal Oxygen "E" Cylinder Holder with IV Pole

STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea
---------	--	------

Wheel Lock Extension

STDS801GT	Wheel Lock Extension	1/ea
-----------	----------------------	------

One Arm Drive

STDS4056	Viper GT 16", 18" and 20"	1/ea
----------	---------------------------	------

8" Caster

STDS1002	Standard 8" Caster	1/ea
----------	--------------------	------

Amputee Adapters

STDS828GT	Amputee Adapters	1/pr
-----------	------------------	------

Flat-Free Wheel with Handrim

STDS840	Flat-free Wheel With Handrims	1/ea
---------	-------------------------------	------

Overhead Anti-Theft Device

STDS823	Overhead Anti-Theft Device	1/ea
STDS834	Overhead Anti-Theft Device Single Pole	1/ea

Wheelchair Cane/Crutch Holder

STDS860	Universal Forward Stabilizer	1/pr
---------	------------------------------	------

Heel loops

STDS831	Heel Loops	1/bx
---------	------------	------

Limb Support

WASR	Limb Support, Right	1/ea
WASL	Limb Support, Left	1/ea

Super Hemi Kit

STDSSHK	Super Hemi Kit	1/ea
---------	----------------	------

Chart Carry Pocket

STDS835	16", 18" and 20" Wheelchairs w/Nylon Upholstery	1/ea
---------	---	------