

Bariatric Sentra EC Heavy-Duty Wheelchair Dual Cross Brace





Item #	Seat Width	Armrest	Product Front Riggings	UOM
STD20ECDDAHD-SF	20"	Detachable Desk Arm	Swing-away Footrests	1/cs
STD20ECDDAHD-ELR	20"	Detachable Desk Arm	Elevating Leg rests	1/cs
STD20ECDFAHD-SF	20"	Detachable Full Arm	Swing-away Footrests	1/cs
STD20ECDFAHD-ELR	20"	Detachable Full Arm	Elevating Leg rests	1/cs
STD22ECDDA-SF	22"	Detachable Desk Arm	Swing-away Footrests	1/cs
STD22ECDDA-ELR	22"	Detachable Desk Arm	Elevating Leg rests	1/cs
STD22ECDFA-SF	22"	Detachable Full Arm	Swing-away Footrests	1/cs
STD22ECDFA-ELR	22"	Detachable Full Arm	Elevating Leg rests	1/cs
STD24ECDDA-SF	24"	Detachable Desk Arm	Swing-away Footrests	1/cs
STD24ECDDA-ELR	24"	Detachable Desk Arm	Elevating Leg rests	1/cs
STD24ECDFA-SF	24"	Detachable Full Arm	Swing-away Footrests	1/cs
STD24ECDFA-ELR	24"	Detachable Full Arm	Elevating Leg rests	1/cs

- Weight Capacity: 450 lbs.
- Carbon Steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish.
- Double embossed vinyl upholstery is durable, lightweight, attractive and easy to clean.
- Standard dual axle provides easy transition of seat to hemi-level.
- Urethane tires mounted on composite wheels provide durability and low maintenance.

- Chrome hand rims standard.
- 8" front casters, adjustable in three positions.
- No maintenance, sealed, precision bearings.
- Padded armrests provide additional comfort.
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings.
- Aluminum foot plates.
- Comes with push-to-lock wheel locks.

28"	Product Weight with 20" Seat	56 lbs.
30"	Product Weight with 22" Seat	57 lbs.
32"	Product Weight with 24" Seat	58 lbs.
12.5"	· · · · · · · · · · · · · · · · · · ·	62 lbs.
	11 0 0	63 lbs.
	5 5	64 lbs.
	11 3 3	
18"	Weight Limit	450 lbs.
42"	Warranty	Limited 3 Year
	30" 32" 12.5" 36" 17.5" - 19.5" 16" 18"	30" Product Weight with 22" Seat 32" Product Weight with 24" Seat 12.5" Carton Shipping Weight with 20" Seat 17.5" - 19.5" Carton Shipping Weight with 22" Seat 18" Weight Limit

accessories

STDS804

Telescoping IV Pole Attachment

STDS820	Telescoping IV Pole Attachment	1/ea
Universal Forw	vard Stabilizer	
STDS862	Universal Forward Stabilizer	1/pr
Anti-Tippers		
STDS818	Anti-Tippers without Wheels	1/pr
Wheel Lock Ex	tension	
STDS801	Wheel Lock Extension	1/ea
Seat Belts		
STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric	1/ea
Universal Oxyg	gen "E" Cylinder Holder	
STDS803	Universal Oxygen "E" Cylinder Holder	1/ea

Universal Oxygen "E" Cylinder Holder with

Overhead Anti-Theft Device

0.0000000000000000000000000000000000000			
STDS823	Overhead Anti-Theft Device		
STDS821	Overhead Anti-Theft Device With IV Hooks		1/ea
STDS834	Overhead Anti-Theft Device Single Pole		1/ea
Anti-Fold Bar			
STDS806	Anti-Fold Bar	20" and 22" Only	1/ea
Wheelchair C	ane/Crutch Holder		
STDS1034	Wheelchair Cane/Crutch Holder		1/ea
Heel Loops			
STDS831	Heel Loops		1/pr
Limb Support			
WASR	Limb Support, Right		1/ea
WASL	Limb Support, Left		1/ea
Chart Carry P	ocket		
STDS836	Chart Carry Pocket	20" Wheelchairs	1/ea

IV Pole Attachment