

## Deluxe Sentra Full Reclining Wheelchair

Dual Axle Wheelchair



HCPCS  
K0001 / E1226

Comes standard  
with swing-away,  
elevating legrests.

- Hydraulic reclining mechanism allows for infinite adjustments up to 180°.
- Durable, heavy-gauge naugahyde upholstery.
- Carbon Steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish.
- Composite Mag-style wheels with chrome hand rims are lightweight and maintenance free.
- 8" front casters.
- Padded armrests provide added patient comfort.
- Wheels set back on frame prevent tipping.
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability.
- Rear anti-tippers standard.
- Comes standard with swing-away elevating legrests.
- Aluminum foot plates are attractive, crack-proof and lightweight.
- Front caster forks are adjustable in three positions.
- Carry pocket standard.
- Headrest extension with cushioned head immobilizer standard.
- Comes with push-to-lock wheel locks.

Item #	Seat Width	Armrest	Product Front Riggings	UOM
STD14RBDDA	14"	Detachable Desk Arm	Swing-away Elevating Leg rests	1/cs
STD16RBDDA	16"	Detachable Desk Arm	Swing-away Elevating Leg rests	1/cs
STD16RBDFDA	16"	Detachable Full Arm	Swing-away Elevating Leg rests	1/cs
STD16RBADDA	16"	Detachable and Adj. Height Desk Arm	Swing-away Elevating Leg rests	1/cs
STD16RBADFA	16"	Detachable and Adj. Height Full Arm	Swing-away Elevating Leg rests	1/cs
STD18RBDDA	18"	Detachable Desk Arm	Swing-away Elevating Leg rests	1/cs
STD18RBDFDA	18"	Detachable Full Arm	Swing-away Elevating Leg rests	1/cs
STD18RBADDA	18"	Detachable and Adj. Height Desk Arm	Swing-away Elevating Leg rests	1/cs
STD18RBADFA	18"	Detachable and Adj. Height Full Arm	Swing-away Elevating Leg rests	1/cs
STD20RBDDA	20"	Detachable Desk Arm	Swing-away Elevating Leg rests	1/cs
STD20RBDFDA	20"	Detachable Full Arm	Swing-away Elevating Leg rests	1/cs
STD20RBADDA	20"	Detachable and Adj. Height Desk Arm	Swing-away Elevating Leg rests	1/cs
STD20RBADFA	20"	Detachable and Adj. Height Full Arm	Swing-away Elevating Leg rests	1/cs
STD22RBDDA	22"	Detachable Desk Arm	Swing-away Elevating Leg rests	1/cs

### specifications

Width Open with 14" Seat	23"	Length with Riggings	51"
Width Open with 16" Seat	25"	Product Weight with 14" Seat	56 lbs.
Width Open with 18" Seat	27"	Product Weight with 16" Seat	59 lbs.
Width Open with 20" Seat	28"	Product Weight with 18" Seat	60 lbs.
Width Open with 22" Seat	30"	Product Weight with 20" Seat	61 lbs.
Width Closed (with 14", 16", 18" 20", 22" seat)	12.5"	Product Weight with 22" Seat	62 lbs.
Height with 14" Seat	48"	Carton Shipping Weight with 14" Seat	65 lbs.
Height with 16" Seat	50"	Carton Shipping Weight with 16" Seat	69 lbs.
Height with 18" Seat	50"	Carton Shipping Weight with 18" Seat	70 lbs.
Height with 20" Seat	50"	Carton Shipping Weight with 20" Seat	71 lbs.
Height with 22" Seat	51"	Carton Shipping Weight with 22" Seat	72 lbs.
Seat-To-Floor Height	17.5" - 19.5"	Weight Limit with 14", 16", 18" and 20" Seat	300 lbs.
Back Height	23"   33"	Weight Limit with 22" Seat	450 lbs.
Seat Depth (with 14", 16", 18", 20" seat)	16"	Warranty	Limited Lifetime
Seat Depth (with 22" seat)	18"		

### accessories

#### Universal Forward Stabilizer

STDS860	Universal Forward Stabilizer	1/pr
---------	------------------------------	------

#### Anti-Tippers

STDS802	Anti-Tippers With Wheels	1/pr
STDS818	Anti-Tippers without Wheels	1/pr

#### Universal O2 Cylinder Holder

STDS6008-1	Oxygen Cylinder Carry Bag for Wheelchair	1/ea
------------	--	------

#### Seat Belts

STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea

#### Oblique Projection Handrim

STDS1022	Oblique Projection Handrim	1/ea
----------	----------------------------	------

#### Standard 8" Caster

STDS1002	Standard 8" Caster	1/ea
----------	--------------------	------

#### Heel Loops

STDS831	Heel Loops	1/bx
---------	------------	------

#### Chart Carry Pocket

STDS835	16", 18" and 20" Wheelchairs w/Nylon Upholstery	1/ea
---------	---	------

#### Limb Support

WASR	Limb Support, Right	1/ea
WASL	Limb Support, Left	1/ea

#### Flat-Free Wheel With Handrims

STDS840	Flat-free Wheel With Handrims	1/ea
---------	-------------------------------	------

#### Light Weight 24" Wheel

STDS1001	Standard 24" Wheel	1/ea
----------	--------------------	------

#### 6" Wheel Lock Extension

STDS801	Wheel Lock Extension	1/ea
---------	----------------------	------

#### Seat Rail Extension

STDSC1818V	Seat Rail Extension Kit	1/pr
STDS811	Seat Rail Extensions	1/pr

#### Overhead Anti-Theft Device

STDS823	Overhead Anti-Theft Device	1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks	1/ea
STDS834	Overhead Anti-Theft Device Single Pole	1/ea

#### Anti-Fold Bar

STDS806	Anti-Fold Bar	1/ea
---------	---------------	------

#### Wheelchair Cane/Crutch Holder

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------