



LC-525i shown in
Aria® Cool Grey fabric

MADE IN
USA
OF FOREIGN & DOMESTIC COMPONENTS

USB
Charger
Port



INFINITY COLLECTION

The Infinity Collection from Pride® Power Lift Recliners features true infinite-position seating and pocketed coils designed to deliver comfort, style and performance. Available in four different fabric selections, the Infinity Collection models deliver maximum relaxation and will match any décor.

AMERICA'S #1
POWER LIFT RECLINER

POWER
RECLINERS

Protect your
investment with
the Pride® Lift
Chair Care Kit



US: 1-800-800-8586 • CAN: 1-888-570-1113
pridemobility.com |

Pride
Mobility Products Corp.
— LIVE YOUR BEST —

INFINITY COLLECTION

AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS



LC-525i shown in
Ultraleather™ Pecan fabric



SPECIFICATIONS

	Petite Wide LC-525iPW ²	Medium LC-525iM ²	Large LC-525iL ²
	Infinite-Position Zero Gravity Trendelenburg Chaise Lounger 5'2" & Below	Infinite-Position Zero Gravity Trendelenburg Chaise Lounger 5'4" - 6'0"	Infinite-Position Zero Gravity Trendelenburg Chaise Lounger 5'9" - 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	19"	20"	21"
Seat Depth	19"	20"	22"
Seat Width	24"	20"	22"
Top of Back to Seat	28"	29"	30"
Back Style	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
Fully Padded Chaise	Yes	Yes	Yes
Heat & Massage	Optional	Optional	Optional
Seating Type	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs
Distance From Wall	19"	20"	21"
Overall Chair Width	38"	34.5"	37"
Chair Weight	130 lbs.	130 lbs.	130 lbs.
Warranty	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*

INFINITY FABRIC SELECTION



DURASOFT



LEXIS STA-KLEEN® sta*kleen. PERFORMANCE FABRIC



ULTRALEATHER™ Uf ULTRALEATHER



ARIA³ CRYPTON. | super look.



Look for the Quick Ship
Indicator for the shortest
lead times available.

BENEFITS

- Premium furniture grade laminate/hardwood frames for added strength and durability
- Pocketed coils deliver individualized support and relaxation
- Industry-first hand control with USB charger port **(not available with Heat and Massage)**
- Standard lumbar pillow
- Standard head & arm covers
- Large, dual pockets
- Easy to remove back
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 9' cord

OPTIONS*

- Deluxe Heat & Massage
- Residential cleaning kit
- Head pillow
- Protective pad
- Lift chair overlay

**Additional lead time may
apply for some options.**

*Additional charges apply for all options.

1. Batteries not included.
2. Available with Deluxe Heat & Massage only. Standard Heat & Massage not available.
3. Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

*Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

US: 1-800-800-8586 • CAN: 1-888-570-1113
pridemobility.com |

Pride
Mobility Products Corp.
— LIVE YOUR BEST —