

Jazzy®

ELITE 14

The Jazzy® Elite 14 is designed for excellent indoor maneuverability and outstanding outdoor capability. The Jazzy Elite 14 features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns around corners. Large front wheels offer extra absorption over obstacles for superb climbing capabilities along with front anti-tip wheels for transitioning over grade changes and handling various terrains with ease.



*Jazzy® Elite 14
shown in Jazzy Red
with high-back seat*

SPECIFICATIONS

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	3"
Turning Radius	24"
Overall Length	35.5"
Overall Width	24.25"
Drive Wheels	14"
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	50A, Dynamic Shark controller
Battery Charger	Off-board, 5 amp
Per-Charge Range	Up to 9.8 miles
Battery Requirements	(2) 12 volt, deep cycle, NF-22
Battery Weight	37 lbs. each
Base Weight	98.5 lbs.
Standard Seat Weight	51.5 lbs. (high-back)
Standard Seat	18" x 18-20" and 20" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

COLOR CHOICES



Jazzy
Red



Jazzy
Blue

OPTIONAL SEATING

- Synergy® Seat
- Solid seat pan

MODELS

Jazzy Elite 14-2S-SS	Code: K0822
Jazzy Elite 14-2S-C	Code: K0823

FEATURES

- Color-through shroud
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick

RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Olympian
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com

Pride[®]
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.