

Jazzy[®]

ELITE ES

The Jazzy[®] Elite ES delivers a reliable blend of power, performance and style for your active lifestyle. The in-line, front-wheel drive technology provides tight turns around corners for optimal indoor maneuverability in small spaces and large front wheels for extra absorption over obstacles. With superb climbing capabilities combined with front anti-tip wheels for transitioning over grade changes easily, the Jazzy Elite ES handles various terrains with ease.



*Jazzy[®] Elite ES
shown in Jazzy Red
with high-back seat*

SPECIFICATIONS

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	0.75"
Turning Radius	23.5"
Overall Length	34.75" without front riggings
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	5" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 12 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	24.5 lbs. each
Base Weight	68.5 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

COLOR CHOICES



Jazzy
Red



Jazzy
Blue

OPTIONAL SEATING

- Solid seat pan

MODELS

Jazzy Elite ES-2S-SS	Code: K0822
Jazzy Elite ES-2S-C	Code: K0823

FEATURES

- Color-through shroud
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap Belt (50", 70", 80", 90")
- Elevating legrest
- Swingaway legrest
- Heel loops
- Military patches
- Height-adjustable footrest bracket
- Swing-away joystick

RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander
- Olympian
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com

Pride[®]
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.