

# Jazzy<sup>®</sup>

## ELITE HD

The Jazzy<sup>®</sup> Elite HD features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front anti-tip wheels for handling various terrains with ease.



*Jazzy<sup>®</sup> Elite HD  
shown in Jazzy Red  
with high-back seat*

## SPECIFICATIONS

<b>Weight Capacity</b>	450 lbs.*
<b>Maximum Speed</b>	Up to 4 mph
<b>Ground Clearance</b>	3"
<b>Turning Radius</b>	24"
<b>Overall Length</b>	35.5"
<b>Overall Width</b>	24.25"
<b>Drive Wheels</b>	14" solid
<b>Front Wheels</b>	3" solid anti-tips
<b>Rear Wheels</b>	6" solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, front-wheel drive
<b>Braking System</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard Electronics</b>	60 amp, Dynamic Shark controller
<b>Battery Charger</b>	Off-board, 5 amp
<b>Per-Charge Range</b>	Up to 9.8 miles
<b>Battery Requirements</b>	(2) 12 volt, deep cycle, NF-22
<b>Battery Weight</b>	37 lbs. each
<b>Base Weight</b>	98.5 lbs.
<b>Standard Seat Weight</b>	51.5 lbs. (high-back)
<b>Standard Seat</b>	20" x 18-20"
<b>Max. Seat Dimensions</b>	24" x 24"
<b>Warranty</b>	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## COLOR CHOICES



Jazzy  
Red



Jazzy  
Blue

## OPTIONAL SEATING

- Synergy® Seat
- Solid seat pan

## MODELS

<b>Jazzy Elite HD-2HD-SS</b>	Code: K0824
<b>Jazzy Elite HD-2HD-C</b>	Code: K0825

## FEATURES

- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 70" lap belt

## OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Safety flag
- Lap belt (50", 60", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Residual limb support
- Heel loops
- Military patches
- Swing-away joystick

## RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Olympian
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

[prideliftsandramps.com](http://prideliftsandramps.com)

**Pride**<sup>®</sup>  
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • [pridemobility.com](http://pridemobility.com)



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.