



LC-358 shown in
Cloud 9 Black Cherry fabric

USB
Charger
Port



MADE IN
USA
OF FOREIGN & DOMESTIC COMPONENTS

HERITAGE COLLECTION

FIRST  **50**
TECHNOLOGY **Faster**
LIFT
SPEED



Scan QR code
for video!

Available on the Heritage Collection Series (except 558, 358XL & XXL)

AMERICA'S #1
POWER LIFT RECLINER

POWER
RECLINERS

Comfort and style complement one another perfectly in the Heritage Collection from Pride® Power Lift Recliners. With an array of plush fabrics to fit any décor and a superior design that provides exceptional strength, Heritage Collection delivers style and durability that will last for years to come. Pride's Heritage Collection offers numerous standard features, superb comfort and quiet and smooth performance, all at a great value.

US: 1-800-800-8586 • CAN: 1-888-570-1113
pridemobility.com | 

Pride
Mobility Products Corp.
— LIVE YOUR BEST —

HERITAGE COLLECTION

AMERICA'S #1
POWER LIFT RECLINER


POWER
LIFT
RECLINERS



	Medium	Petite Wide	Small	Medium	Large	Heavy-duty	Very Heavy-duty
SPECIFICATIONS	LC-558 ⁸	LC-358PW ^{2,5}	LC-358S ^{2,5}	LC-358M ^{2,5}	LC-358L ^{2,5}	LC-358XL ^{4,7}	LC-358XXL ^{4,7}
	Infinite, Lay Flat Chaise Lounger 5'4" - 6'	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4" - 6'	3-Position Chaise Lounger 5'8" - 6'2"	3-Position Chaise Lounger 5'6" - 6'1"	2-Position Chaise Lounger 5'6" - 6'1"
Weight Capacity	325 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	19"	19"	19"	20"	21"	20"	20"
Seat Depth	20"	19"	19"	20"	22"	22"	22"
Seat Width	20"	24"	18"	20"	22"	26"	30"
Top of Back to Seat	26"	23.5"	24"	24.5"	25.5"	29"	29"
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage	No	Optional	Optional	Optional	Optional	Optional	No
Seating Type	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Sinuous Spring	Sinuous Spring
Distance From Wall	24"	18"	18"	18"	19"	21"	14"
Overall Chair Width	33"	39"	35"	35"	37"	39.5"	43"
Chair Weight	147 lbs.	118 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
Warranty	7-Year Prorated ¹	7-Year Prorated ¹	7-Year Prorated ¹	7-Year Prorated ¹	7-Year Prorated ¹	7-Year Prorated	7-Year Prorated

HERITAGE FABRIC SELECTION


QS CLOUD 9




LEXIS STA-KLEEN® sta-kleen® PERFORMANCE FABRIC



ULTRALEATHER™ Uf upholstery



ARIA³ CRYPTON® | super fabric



Look for the Quick Ship Indicator for the shortest lead times available.



1. Batteries not included.
2. Available with footrest extension.
3. Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.
4. Removable, fully padded chaise.
5. Unique spring design and a one-piece soft foam seat for outstanding relaxation.
6. Features Dual Motor Sync System.
7. Only available in Stone, Walnut, Lexis Sta-Kleen, Chestnut and Mushroom.
8. Standard with footrest extension and lumbar pillow.

BENEFITS

- Premium furniture grade laminate/hardwood frame for added strength and durability
- Industry-first hand control with USB charger port (**not available with Heat and Massage**)
- Standard head & arm covers, large dual pockets
- Removable back
- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation (except LC-358XL & LC-358XXL)
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 9' cord

OPTIONS*

- Heat & Massage (except XXL and LC-558)
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Non-skid leg levelers
- Footrest extension available (except LC-358S, LC-358XL and LC-358XXL)
- First Up technology (except LC-358XL, LC-358XXL, and LC-558)
- Lift chair overlay

Additional lead time may apply for some options.

*Additional charges apply for all options.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.
*Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

US: 1-800-800-8586 • CAN: 1-888-570-1113
pridemobility.com |

Pride
Mobility Products Corp.
— LIVE YOUR BEST —